

NLIS CATTLE ORDER FORM

Fax 07 3245 9110 Email: customerservice@allflex.com.au



Retailer Details:

Retailer:

GDL – Blackall
117 Shamrock St
BLACKALL 4472

Contact Person: **Ann Russell**

Email Address: **arussell@gdlrural.com.au**

Purchase Order Number:

Deliver to:


Retailer Applicant

Applicant Details:

07 4657 6000

Applicant Name						Date		
Trading Name				Property Identification Code				
Property Address								
Signature				Phone Number				
				Email Address				

Tag Order Details:

Description	Series Start Number	Quantity
White Breeder Tag: Cattle born on property. If you have previously ordered Allflex NLIS devices your serial numbers will automatically follow on from the last order. If you require a different series of numbers enter the new start number.		
Orange Post-Breeder Tag: Cattle brought onto the property or Replacement Tag. If you have previously ordered Allflex NLIS devices your serial numbers will automatically follow on from the last order. If you require a different series of numbers enter the new start number.		
Matching Management Tags: Includes male button. Marking is PIC over Series Number. Note: Tags not shown to scale Colour: <input type="checkbox"/> Yellow <input type="checkbox"/> Pink <input type="checkbox"/> Green <input type="checkbox"/> White <input type="checkbox"/> Orange <input type="checkbox"/> Blue <input type="checkbox"/> Red <input type="checkbox"/> Purple <input type="checkbox"/> Black  Maxi: <input type="checkbox"/> Quantity: <input type="text"/> Series Start Number: <input type="text"/> LARGE: <input type="checkbox"/> Quantity: <input type="text"/> Series Start Number: <input type="text"/> Quantity: <input type="text"/>		
Customised Matching Management Tags Includes Male Button. Available in the same colours as above. Any marking variation outside of PIC over Series Number is a Customised Tag. Please indicate the colour and size required (Large or Maxi) and tag marking layout above.		
DNA Collector Pre-printed DNA sample collector linked to NLIS device and/or matching visual tags. Self adhesive collector for hair root samples to be collected at time of tag application and archived.		
Applicator Universal Applicator (red) or IdentiPlier (green) are used for both NLIS devices and Matching Visual Tags. Alternatively the UltraMatic Applicator can be used for NLIS devices or LazaMatic Applicator can be used for Visual Tags.		
Freight Standard Australia Post (Default) <input type="checkbox"/> TNT <input type="checkbox"/> Express Post <input type="checkbox"/> Other... <input type="text"/>		